

## LABORATORY HEATING MANTLE



- Rugged, steel exterior.
- Heating elements are flexible and elastic ensuring close contact with the flask surface minimizing breakage of glassware.
- Glasswool insulation ensures minimal heat loss.
- Uniform surface temperature upto 450°C.
- Electrically operated on 230 VAC, Single Phase.

Model Code	Flask Size	Rating
HHM-3	250ml	150 W
HHM-4	500ml	200 W
HHM-5	1 Lt.	300 W
HHM-6	2 Lt.	450 W
HHM-7	3 Lt.	450 W
HHM-8	5 Lt.	600 W

**HARIDAS INSTRUMENTS & EQUIPMENTS CO.**

**(An ISO 9001:2008 Company)**

Sector-23, Plot No.76, CIDCO Service Industrial Area, Turbhe, Navi Mumbai-400 705, Maharashtra, India.  
 Tel : 91-022-27833046, 27835042 Fax : 91-022-27835042 email : kph@vsnl.com url : [http:// www.haridasinstruments.com](http://www.haridasinstruments.com)