

LABORATORY HOT AIR OVEN



- Rugged, steel construction.
- Interior available with Mild Steel or Stainless Steel.
- Available with or without forced air circulation.
- Temperature range upto 250°C.
- Choice of thermostatic , digital On/Off or PID temperature control.
- Electrically operated on 230 VAC, Single Phase.

Model Code	Chamber Size (H x W x D)	Shelves	Rating
HLO - 30	30x30x30cm	1	700 W
HLO - 35	35x35x35cm	2	750 W
HLO - 45	45x45x45cm	2	1.2 kW
HLO - 54	60x45x45cm	3	1.7 kW
HLO - 60	60x60x60cm	3	2 kW
HLO - 65	90x60x45cm	3	2.5 kW
HLO - 90	90x60x60cm	3	2.5 kW

HARIDAS INSTRUMENTS & EQUIPMENTS CO.
(An ISO 9001:2008 Company)

Sector-23, Plot No.76, CIDCO Service Industrial Area, Turbhe, Navi Mumbai-400 705, Maharashtra, India.
 Tel : 91-022-27833046, 27835042 Fax : 91-022-27835042 email : kph@vsnl.com url: http:// www.haridasinstruments.com